## MADiE 1.1.14 Release Notes

The following updates are included in the MADiE 1.1.14 release:

Category	Feature Description
VSAC Authentication	Replaced the authentication method used when querying VSAC. The TGT/ST workflow will be replaced with Basic Authentication header containing the user's VSAC API Key.
Measure Creation	Enhanced CQL Editor error messaging so that warnings can be distinguished from errors and warnings do not block test case execution. (BONNIEMAT-1353)
Measure Export	<ul> <li>Updates have been applied to Measure exports so that they are formatted correctly. The updates include the following: <ul> <li>The fullURL for all library type resources to have a root of "<u>https://madie.cms.gov</u>" instead of "<u>http://ecqi.healthit.gov</u>"</li> <li>Set the measure type resource URL element to the same value as the fullURL value (Measure.url = Measure.fullURL)</li> <li>Changed the PUT requests on all resources to POST requests (change "method": "POST" to "method": "PUT")</li> <li>In the Library resources (primary and all included libraries) updated the relatedArtifact array: <ul> <li>Replace "url" with "resource"</li> <li>Include the version of the resource when applicable at the end of the resource value (pipe-delimited) for measure model and libraries, NOT for value sets or code systems</li> </ul> </li> </ul></li></ul>
Test Case Import	<b>QI Core Measures:</b> Warning message has been updated so that when a test case that does not have a unique test case name (the family and given name combination across the measure) is imported, the message will indicate exactly what the application considers non-unique. (BONNIEMAT-1602)
Test Case Export	<b>QI Core Measures:</b> MADiE has been updated to add a 'request' field when the test case is Transaction bundle type. When a user sets the test case bundle type as "Transaction", the export will now include the following:
	<pre>"method": "POST", "url": "MeasureReport/<measure id="" report="">" }</measure></pre>
-	(BONNIEMAT-1607)
Test Case Creation	<b>QI Core Measures:</b> MADiE has been updated so that the reference field(s), for patient references, will automatically be set to the format of "Patient/ <patientid>", when a user saves a test case. For example, reference fields that were in the format "Patient/<anything>" will be saved in the format of "Patient/<patientid>" (BONNIEMAT-1608)</patientid></anything></patientid>