

MADiE 2.0.2 Release Notes

The following updates are included in the MADiE 2.0.2 release:

Category	Feature Description
Test Cases Creation	MADiE has been updated to allow users to shift test case dates either forwards or backwards, a user-specified number of years. This functionality will be available for a single test case or for all test cases on a measure.
Measure Export	<p>QI-Core Measures: The way MADiE passes values to the liquid template has been updated to utilize the display name instead of the code name for identifier value sets, so that the values will appear in sentence-style case and be easier to read. (MADIE-1729)</p> <p>Measures transferred from MAT to MADiE without a Stratification Description now displays "None" in the Stratification Description in the human readable export. (MADIE-1959)</p>
Measure Creation	<p>Users can now associate a CMS ID from a QDM measure to a QI-Core measure, indicating they are the same measure. Users can also copy over the metadata from the QDM measure to the QI-Core measure. (MADIE-1900)</p> <p>Users now have the ability to request CMS IDs be deleted from their draft measure when there are no prior versions of the measure. To request that a CMS ID be deleted from your measure, fill out the MADiE CMS ID Deletion Request Form and send it to semanticbits-madie-help@icf.com.</p> <p>When a user clicks the Generate CMS ID link on the Measure Details screen, a popup will display asking users to confirm they want to generate a CMS ID for the measure.</p> <p>Users can now request to have an unused library they own to be deleted after verifying no other measures or libraries include that library by submitting that request to the MADiE Issue Tracker. All versions (including draft) of the library will be deleted. (MADIE-1801)</p> <p>QDM Measures: MADiE has been updated to allow users to include 128 SDEs and RAVs in a measure. (MADIE-1921)</p> <p>QI-Core Measures: MADiE now requires that all QI-Core measures and libraries include the FHIRHelpers library. If you do not include the FHIRHelpers library in your measure or library, you will receive a CQL error.</p>
UI Enhancements	Users can now move any popup screen around to see the contents of the screen behind it. (MADIE-1965)
QMIG STU4 Updates	Users can now select values for "Include in Report Type" for each SDE and RAV on their measures.

MADiE 2.0.2 Release Notes

Category	Feature Description
Resolved Known Issues	<p>When users are adding SDEs and RAVs users some of the available options listed in the drop down can't be viewed. (MADIE-1969)</p> <p>Date Added: 12/21/2022</p> <p>When users search for measure names and include special characters in the search criteria, the user may receive a 403 error, or the search results will not be returned.</p> <p>Date Added: 6/26/2024</p> <p>QDM Measures: In some cases, large value sets were not being fully retrieved. As a result, test cases may fail unexpectedly as some codes may not be pulled in. (MADIE-1950)</p> <p>Date Added: 6/26/2024</p> <p>QDM Measures: Edge test cases may fail unexpectedly due to the measurement period end date precision being used in MADiE. (MADIE-1953)</p>