

## MADiE Known Issues

### CQL-to-ELM Translator Known Issues

Date Added	Issue Description
4/12/2023	<p>There is an issue in how reasonCode is getting mapped to ELM for the CommunicationNotDone datatype in QI-Core 4.1.1. Refer to the following GitHub ticket for more information on this issue:  <a href="https://github.com/cqframework/clinical_quality_language/issues/1153">https://github.com/cqframework/clinical_quality_language/issues/1153</a>  <a href="#">BONNIEMAT-1420</a></p>

### MADiE Known Issues

Date Added	Issue Description
10/10/2023	<p>Users will automatically be logged out of MADiE after 30 minutes of inactivity. If the user has more than one active session of MADiE in any combination (in different tabs in the same browser or in different browsers) and one of the sessions reaches 30 minutes of inactivity the session will expire, and the user will be logged out of all sessions. If the user is logged out prior to saving changes, any unsaved changes will not be retained when the user logs into MADiE again.</p> <p><b>Work around:</b> Users can keep all sessions active by clicking on the screen for each session periodically. Clicking on each active session prior to reaching 30 minutes of inactivity will prevent the sessions from expiring and prevent the user from being logged out of all sessions.</p> <p><a href="#">BONNIEMAT-1554</a>  <a href="#">BONNIEMAT-1623</a></p>
11/15/2023	<p><b>QDM Measures:</b> Status' date timing included in QDM test case JSON is not being saved after translation completes.</p>
3/27/2024	<p>Measures transferred from MAT to MADiE need to be reviewed in MADiE and the CQL updated (can just add a space to the CQL) after the transfer and prior to working with the measure (including working with test cases or exporting the measure). Until this step is complete you may see a CQL error notification on the Test Cases screen and/or during Measure Export.</p>
4/19/2024	<p><b>QDM Measures:</b> When editing a MADiE test case that was imported from Bonnie QDM, if the test case description exceeds the allowable 250 characters, the 'Save' button is disabled and updates cannot be made to the test case.</p> <p><b>Work around:</b> Edit the MADiE test case description to less than 250 characters.</p> <p><a href="#">MADIE-1825</a></p>

<b>Date Added</b>	<b>Issue Description</b>
6/12/2024	<p>Users may experience toggling back and forth between their current tab and a new tab or flickering and, in some cases, discontinuation of their MADiE session.</p> <p><b>Work around:</b> Refresh the page or log out of all MADiE instances and begin a new MADiE session.  <a href="#">MADIE-1922</a></p>
8/15/2024	<p>If a MADiE user's UMLS API key has changed, the original UMLS API key entered into MADiE cannot be updated by the user.</p> <p><b>Work around:</b> Open a MADiE Issue Tracker ticket to make a request for the MADiE team to manually remove the existing UMLS key, so a new UMLS key may be entered.  <a href="#">MADIE-1912</a></p>
8/30/2024	<p>When importing QI-Core test cases from MADiE into a measure with Stratifications, the test cases will appear as if the measure does not have stratifications. This occurs even if the measure the test cases were exported from had stratifications.</p> <p><b>Work around:</b> Navigate to the population criteria tab, modify any stratification, click save. Then navigate back to the test cases and stratification should now be present.</p>
10/28/2024	<p><b>QDM Measures:</b> Measures making comparisons on quantities may receive inaccurate test case results. <a href="#">MADIE-2056</a></p>